



Gly/TCA/nucleotide and NAD related metabolites

Service Code: NAD

Summary: Profile of Central Metabolism, including glycolysis, pentose-phosphate shunt, TCA cycle, nucleotide pool and NAD related metabolites. One step organic solvent extraction of cultured cells or tissues, separated on a 1mm x150mm HILIC column in a 35 min cycle. All analytes and Internal Standards are measured by ESI⁻ ionization on a LC-QTOF mass spectrometer and reported as uM normalized to wet tissue weight or cell proteins. CV's are generally 15%.

Container: Eppendorf Tube or equivalent

Normal Volume: Plasma (100 ul) Tissue (50-100 mgs); Cells (2E7).

Minimal Volume: Plasma (50 uL)Tissue (30 mg); Cells (~5E6)

Special Handling: If human or primate, note any known presence of infectious agents.

Sample Collection: Snap freeze by liquid nitrogen. For tissues, resect and snap-freeze as soon as practical in tared centrifuge tube. Provide both sample weight and tared vial weight on sample submission

Reference: [Matthew A. Lorenz](#), [Charles F. Burant](#), and [Robert T. Kennedy](#) (2011) "Reducing Time and Increasing Sensitivity in Sample Preparation for Adherent Mammalian Cell Metabolomics", *Anal. Chem.* 83(9): 3406–3414.

Table I: Analytes reported:

Analyte	Abbr.	Mol Formula	LOQ(uM)
Acetyl-CoA	aCoA	C ₂₃ H ₃₈ N ₇ O ₁₇ P ₃ S	1
Citrate/Isocitrate combined	Cit/i-Cit	C ₆ H ₈ O ₇	1
Succinate	Suc	C ₄ H ₆ O ₄	1
Malate	Mal	C ₄ H ₆ O ₅	1
2-Phosphoglycerate/3-Phosphoglycerate combined	2PG/3PG	C ₃ H ₇ O ₇ P	1
Phosphoenolpyruvate	PEP	C ₃ H ₅ O ₆ P	1
Adenosine monophosphate	AMP	C ₁₀ H ₁₄ N ₅ O ₇ P	1
Adenosine diphosphate	ADP	C ₁₅ H ₂₃ N ₅ O ₁₄ P ₂	1
Adenosine triphosphate	ATP	C ₁₀ H ₁₆ N ₅ O ₁₃ P ₃	1
Flavin adenine dinucleotide	FAD	C ₂₇ H ₃₃ N ₉ O ₁₅ P ₂	1
Nicotinamide adenine dinucleotide	NAD	C ₂₁ H ₂₈ N ₇ O ₁₄ P ₂	1
Nicotinamide adenine dinucleotide, reduced	NADH	C ₂₁ H ₂₉ N ₇ O ₁₄ P ₂	1
Nicotinamide adenine dinucleotide phosphate	NADP	C ₂₁ H ₂₉ N ₇ O ₁₇ P ₃	1
Nicotinamide adenine dinucleotide phosphate, reduced	NADPH	C ₂₁ H ₃₀ N ₇ O ₁₇ P ₃	1



MEDICAL SCHOOL
METABOLOMICS CORE
UNIVERSITY OF MICHIGAN

Erythrose 4-phosphate*	E4P	C ₄ H ₉ O ₇ P	1
Ribulose 5-phosphate/Xylulose 5-phosphate/ribose-5-phosphate combined*	R5P/X5P/ Ru5P	C ₅ H ₁₁ O ₈ P	1
6-phosphogluconate*	6PG	C ₆ H ₁₃ O ₁₀ P	1
Sedoheptulose 7-phosphate*	S7P	C ₇ H ₁₅ O ₁₀ P	1
Fructose-6-phosphate + glucose-6-phosphate	F6P/G6P	C ₆ H ₁₃ O ₉ P	1
Fructose-bisphosphate	FBP	C ₆ H ₁₄ O ₁₂ P ₂	1
Nicotinic acid (NA)	NA	C ₆ NH ₅ O ₂	1
Nicotinic acid mononucleotide (NaMN)	NaMN	C ₁₁ H ₁₅ N ₂ O ₈ P	1
Quinolinic acid (QA)	QA	C ₇ H ₅ NO ₄	1
nicotinamide mononucleotide (NMN)	NMN	C ₁₁ H ₁₅ N ₂ O ₈ P	1
nicotinamide (NAM)	NAM	C ₆ H ₆ N ₂ O	1

*Metabolites are low concentrations and below detection limit in some samples.